

Hela Cleaning-Kit 9050

Hela Cleaning Gun 9011

Stainless Steel

Deionized Water

Sterilizable

The Stainless Steel Cleaning-Kit comes with 6 screwable attachments. The set was specially developed for the use of demineralized water. The Cleaning Gun is easy to take apart by hand without using tools, this allows an effortless cleaning of the cleaning gun. The use of high-quality materials in combination with latest manufacturing engineering ensure an unrestricted steam sterilizability. The cleaning gun also makes an ergonomic work possible in use of high working pressure in a wide temperature range.

Premium Set 9052

- ① 1115 Hose attachment
- ② 1040 Pointy attachment for the cleaning of valves
- ③ 1120 Luer attachment for syringes and cannula with Luer connector
- ④ 1125 Record attachment for syringes and cannula with record cone
- ⑤ 1130 Attachment for the cleaning of drainage tubes
- ⑥ 5050 Spray jet for the cleaning of glassware and surfaces



Attachment holder made out of stainless steel for the attachment storage at the sink or for a save wall mounting.