




# Make it visible!

PEREG Cleaning Efficacy Indicators for validation and routine monitoring according to EN ISO 15883



## TOSI®

Monitoring the cleaning efficacy of washer disinfectors.  
Ready made test object with test soil correlating to blood.



## TOSI®-LumCheck System

Monitoring the cleaning efficacy of cannulated instruments.



## TOSI®-FlexiCheck System

Monitoring the cleaning efficacy of flexible endoscope washers.



## SonoCheck

Monitoring the energy of ultrasonic cleaner.



## HemoCheck-S

Test kit for detection of blood-protein residue.  
A sensitive test method with a detection limit of 0.1µg.



## Pyromol-E

Test kit for detection of protein residue in flexible endoscopes.



## HydroCheck-E

*A new innovation by PEREG*

A test for residual water in flexible endoscopes.